

Capno-T Specifications

Palm Carbon Dioxide Monitor
Mainstream Capnograph





Features

- Small and super lightweight unit with ultimate portability, ideal for emergency, pre-hospital first aid, patient transit, operating room, ICU, and sleep study.
- Large reading display, clearly visible from more than 5m away.
- Turn on the unit and it is ready to use instantly, with no need to calibrate.
- Minimal warm-up time with first Etco2 value displayed in 5 seconds, and accurate readings displayed in 15 seconds.
- Advanced automatic atmosphere pressure compensation, and the gas compensation programs ensures accurate measurement results.
- Data output via USB and Bluetooth (optional) connections, compatible with multiple database software systems and Bluetooth printers.
- Sturdy design, suitable for use in various challenging clinical environments.



Technical specifications

● Power

AC Input:100V~250V, 50Hz/60 Hz

AC Power Consumption: 5W

● Operating environment

Working temperature: 5~40℃

Relative humidity: 30-75% (non-condensing)

Atmospheric pressure: 86~106kpa

Co2/Etco2 Measuring Range: 0~150mmHg

● Accuracy

0~40mmHg, Inaccuracy≤± 2mmHg;

41~70mmHg, Inaccuracy≤±5% of Reading;

70~100mmHg, Inaccuracy≤± 8% of Reading;

101~150mmHg, Inaccuracy≤±10% of Reading;

Respiration Rate: 3~150t/min

Inaccuracy: The bigger value of ±1% or ±1t/min

Battery: Lithium battery, rechargeable. 3.7Volt,1.2AH

Dimension: 38(W)*42(H)*44(D)

Net Weight: 80g

Display size: 1.5 inches



Accessories



A3 Adult Airway adapter



A3N Pediatric Airway adapter

● Applicable Departments

Anesthesiology, operating room, respiratory department, first aid, emergency department, ICU, NICU, and other departments

● Application

Airway positioning, breathing monitoring, circulation monitoring, patient evaluation and patient screening.